

Multi-Parameter Indicator Type 4 VP-Code/Quality: 645_50



1. Product description

Mult-Parameter strips objectively test your sterilisation process according to sufficient steam penetration

2. Sterilisation suitability

Steam sterilisation

3. Normative requirements

The indicator strips meet the requirement of DIN EN ISO 11140 Part 1

4. Statutory requirements

The indicator strips complies with the following regulations and directives:

- Regulation (EG) Nr. 1907/2006
- Directive 2011/65/EU Art. 4 und Annex II
- Directive 94/62/EC
The sums of the concentrations of lead, cadmium, mercury and hexavalent chromium (Cr V) are below 100 ppm.

5. Substances

Our product does not contain natural latex, PVC, Bisphenol A, Colophony, TSE/BSE risk materials, SVHC Substances and Phthalates as part of the formulations.

6. Product properties

Indicator:

Indicator change from blue to pink

Technical Data Sheet



Material:

Inspection characteristic	Typical Value	Unit	Standards
Grammage material	240	g/m ²	DIN EN ISO 536
Lamination OPP	12	µm	

Lamination:

one-sided on printed side

Print:

2-coloured

7. Identification

Carton label



8. Disposal

Indicator strips have to be disposed of as non-recyclable waste (waste code 15 01 06).

9. Assortment

Indicator strips Type 4 (STEAM)

<https://www.stericlin.de>

10. Field of application

This Technical Data Sheet refers to the above-mentioned product group/quality and applies until the next revised edition. Other product-related documents are available on request if required.

11. Additional information

The information contained in this Technical Data Sheet (TDS) is based on our current knowledge and experience.

VP does not guarantee and makes no representation as to the accuracy or completeness of the information contained herein.

VP also assumes no responsibility regarding the consequences arising from the use of information or with regard to any misprints.

The user is responsible for thoroughly checking our products for their suitability for the specific application requirements of the user as well as to their suitability for the respective conditions to which the product is exposed.

The user also has the responsibility for the proper, safe and lawful use, processing and handling of our products, especially when recommendations for safe use and storage have been given. Under no circumstances shall the information contained herein be construed as representation, warranty or guarantee.