



VWR® Pestle for Cell Strainer


Description: Pestle for Cell Strainer, Green, Sterile

Purpose: Supporting the use of cell strainer, reduce sample loss

Materials


Pestle: PP (Polypropylene) Color: Green

Dimensions (Unit: mm)



Features

- Made from high grade gamma resistance polypropylene
- Convex Pestle head includes molded textured surface to better manipulate material in the cell strainer
- A mesh bottom line, grinding better
- Handle special process design, non-slip and easy grip
- DNase/RNase free and non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® Pestle for Cell Strainer

VWR EU EAN. No.	Description	Sterile	Packout Description
734-2819	Pestle for Cell Strainer, Green, Sterile, Individually Packaged	Yes	Individually Packaged, 100/case
734-2820	Pestle for Cell Strainer, Green, Sterile, Bulk	Yes	25/bag, 50/case