

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Rectangular or Square Bottles, Wide Neck with Screw Cap
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2022-02-23

Catalogue Number	Capacity [mL]	Thread	WxDxH [mm]	Type
215-1904	4000	100 - 415	147 × 147 × 290	Square
215-1905	125	28 - 415	38 × 59 × 97	Rectangular
215-1906	250	38 - 415	51.5 × 74.5 × 112	Rectangular
215-1907	500	48 - 415	64 × 99.5 × 114	Rectangular
215-1908	1000	53 - 415	67 × 127 × 180	Rectangular
215-1909	2000	63 - 415	96 × 150 × 237	Rectangular

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-04	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	NA
Dimensional	QSF-QCIR-02	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material	Carboy Closure	High Density Polyethylene (HDPE). Polypropylene (PP).
		Above material is 100% Virgin, Compliance with FDA 21 CFR
Temperature Resistance		up to +120 °C
Sterilization		Not applicable.
ATP Assay		Not applicable.
DNase & RNase Free		Not applicable.
BSE/TSE		No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity		Not applicable.
Non-Pyrogenic Statement		Not applicable.
Latex Statement		Not applicable.
BPA Statement		Not applicable.
DEHP Statement		Please see RoHS information below.
RoHS		No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement		No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.

Storage Conditions Store in Room temperature.

Shelf Life Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature