

Material	20413.420
Material description	Nitric acid 68%
Grade	GPR RECTAPUR
CAS Number	7697-37-2
Molecular formula	HNO ₃
Molecular mass	63.01

Characteristics	Specifications
Assay	68.0 - 70.0 %
Density (20/4)	1.400 - 1.420
Heavy metals (as Pb)	≤ 20 ppm
Ignition residue (SO ₄)	≤ 50 ppm
Cl (Chloride)	≤ 5 ppm
SO ₄ (Sulphate)	≤ 5 ppm
Fe (Iron)	≤ 2 ppm
Pb (Lead)	≤ 1 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.