

Cole-Parmer®

Cole-Parmer® Micro-Kjeldahl Extraction Equipment - EH-200

- Rugged, easy to clean stainless steel construction
- Back-mounted brackets hold flask support rods (supplied with each unit)
- 6 positions which can accommodate flasks of between 18ml - 50ml capacity
- Dedicated controllers for each recess are housed in a “cool zone”
- Grounded stainless steel earth screen covering the heating element



Cole-Parmer® Micro-Kjeldahl Extraction Equipment - EH-200

Cole-Parmer® EH-200 Micro-Kjeldahl equipment has been designed for Kjeldahl extraction on a micro scale. These multibank units have 6 recesses in a stainless steel outer casing with back-mounted brackets (supplied with each unit), which can hold support rods of 12.7mm diameter max. This 6 recess model accommodates flasks of between 18-50mL only.

Each heating position has its own energy regulator incorporating an On/Off switch and a “Mains to Heater” amber neon indicator. There is also a “Mains On” clear neon indicator on the front panel. The lower part of the unit houses dedicated controllers for each recess. This “cool zone” housing is separated from the heating element by a stainless steel screen and a well-ventilated air space. The heating element consists of thermally insulated element wire stitched into a cartridge, and operates in the temperature range of 550°C to 800°C max.

The EH-200 Micro-Kjeldahl equipment models incorporate an earth screen and are double-fused for extra safety. The rugged stainless steel outer casing is durable and easy to clean.



Technical Specification

Specification	EH-200
Heating element temp	550°C - 800°C maximum
Case material	Stainless steel
Thermal insulation	Ceramic fibre/mineral wool
Clamps for support rods	Fitted with adjustable clamps to accept the standard arms supplied
Wells capacity, ml	18-50ml

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
EH-200-6 Micro Extraction Heater	03074-32	EH-200	EH-200-M-6	MM2313/E
EH-200-6-115 Micro Extraction Heater	03074-30	EH-200	EH-200-M-6-115	MM2313/EX1
EH-200-6-EU Micro Extraction Heater	03074-34	EH-200	EH-200-M-6-EU	MM2313/EX6

Accessory Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Fume collection tube for MM series	03074-99	-	-	AT4001

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer®
essentials
an antylia scientific company

Find out more at antylia.com