



ACROS ORGANICS part of Thermo Fisher Scientific

ThermoFisher
S C I E N T I F I C

Product Specifications

Product 416780010
Paraformaldehyde, 96%, extra pure

General Product Data

Version	4
CAS No	30525-89-4
Molecular weight	0.00
Molecular formula	
Linear formula	(CH ₂ O) _n
Flash point (°C)	71

Product Specifications

Appearance (Color)	White
Appearance (Form)	Powder or beads
Infrared spectrum	Conforms
Titration with HCl	>=95 %
Ash	=<0.02 %
pH	4.0 to 5.5 (10 % suspension)
Free acid (HCOOH)	=<0.03 %

Product Specifications are subject to amendment and may change over time.

Acros Organics BV

ENA23, zone 1, nr 1350, Janssen Pharmaceuticaaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: <http://www.thermofisher.com>
1 Reagent Lane, Fair Lawn, NJ 07410, USA
Fax 201-796-1329