

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Spreaders, Triangle Shaped, Disposable
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2025-12-05

Catalogue Number	Sterility	Colour	Packed	PK
612-2688	Sterile	Blue	Individually wrapped	25
612-2690	Sterile	Yellow	Individually wrapped	25
612-2687	Non-sterile	Blue	Bulk	25
612-2689	Non-sterile	Yellow	Bulk	25

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared the products are manufactured in accordance with their internal quality system. Products are inspected and controlled through whole production processing in accordance with their applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene.
Sterilization	Products are E-beam irradiated with a specified dose range of 25-25 kGy. (Which is not applicable for non-sterile products).
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE) according to EMEA 410/01 revision 3.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Product is pyrogen free.
Latex Statement	Product is free from latex.
BPA Statement	Bisphenols are not present.
DEHP Statement	DEHP not present.
RoHS	Product conforms to the requirements of the European Union's Restriction On Use Of Certain Hazardous Substances In Electrical And Electronic Equipment (RoHS) Directives 2011/65/EU and 2015/863/EU (RoHS3).
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.
Storage Conditions	Dry, Room Temperature.
Shelf Life	36 Months.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	12-2025