

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Combined pH Electrodes, pHASE
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2026-03-04

0-14 pH, BNC Connection, 12 mm diameter, 150 mm length

Catalogue Number	Type	Material	Electrolyte	For
662-1782	Refillable, double junction (plastic)	Epoxy	3.5 mol/l KCl	Most applications
662-1800	Refillable, double junction (plastic)	Glass	3.5 mol/l KCl	Aqueous media, beverage, Tris buffer, cosmetics, ultra-pure water
662-1778	Sealed, double junction (Pellon)	Epoxy	Gel	Aqueous media, beverage, Tris buffer
662-1799	Sealed, double junction (Pellon)	Glass	Gel	Aqueous media, beverage, Tris buffer
662-1780	Sealed, flat tip, double junction (HDPE)	Epoxy	Gel	Aqueous media, beverage, Tris buffer, cosmetics, ultra-pure water

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative product sensors are individually tested according to quality procedure DOP-Q126.

Product Specifications

Material Glass, Ceramic.

Temperature range 0 to +100 °C

VWR International bv, Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium - <http://www.vwr.com>

Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not Applicable.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Product contains DEHP.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Product meets the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006.
Storage Conditions	Store in supplied bottle of solution in cool dry place.
Shelf Life	1 year from manufacture date.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	03-2026