



CERTIFICATE OF CONFORMANCE AND ANALYSIS

Date: 08/01/2019

This certification is provided as full assurance that the following product code was manufactured in accordance with prescribed procedures and specifications.

Product Code:	RE-1284F
Description:	Sporicidin Disinfectant Solution
Lot Number:	3275129C
Batch Size:	2970 Gallons
US EPA Registration Number:	8383-3
Manufacturing Date:	02/26/2019
Expiration Date:	02/2023

Description	Specification	Method	Results
Appearance	Match Std	Visual	Pass
Odor	Match Std	Olfactory	Pass
pH (AS IS)	8.3-8.7	pH Meter/ AS IS	8.5
Total Phenol	1.54 – 1.70%	HPLC or Titration	1.61%
Specific Gravity	Record	Pycnometer or Hydrometer	1.007

Sample Conforms: YES

This product meets the requirements of Consumer Product Safety Improvements Act of 2008

A handwritten signature in black ink, appearing to read "C. Inthavongsa".

Champagne Inthavongsa
Quality Supervisor
cinthavongsa@contecinc.com