

**Technical Data Sheet Nr. 3871**

Product Code : 3871. Mojonnier extraction flask for Fat Analysis

Batch-No: LIE320158

Applications: It is used to measure the fat content especially in dairy products. The methodology is described in ISO 1211, ISO 1735, ISO 1736, ISO 1737; ISO 1854; ISO 2450, ISO 5543.

Features: The tube is made from clear borosilicate glass 1,4mm in thickness. The precise dimensions of the tube comply with the dimensional requirements of ISO 3889 as described in Table 1. The neck of the flask allows the tight closure with a stopper

International Standard Compliance: The product complies with ISO 3889