

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Spinbar® Magnetic Stir Bars, Egg-Shaped
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2024-05-17

VWR Catalogue Number	Length in mm [inches]	Diameter in mm [inches]
58949-012	9.5 mm (3/8")	5 mm (3/16")
58949-010	15 mm (5/8")	6 mm (1/4")
58949-006	19 mm (3/4")	9.5 mm (3/8 ")
58949-196	25 mm (1")	12 mm (1/2")
58949-016	38.1 mm (1 1/2")	15.9 mm (5/8 ")
58949-200	32 mm (1 1/4")	15.9 mm (5/8 ")
58949-210	41 mm (1 5/8")	19 mm (3/4 ")
58949-002	51 mm (2")	19 mm (3/4 ")
58949-243	64 mm (2 1/2")	19 mm (3/4 ")
58949-008	76 mm (3")	19 mm (3/4 ")

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Autoclavable at 121° C for 20 minutes.

Tolerance ± 0.020"

Product Specifications

Material	Dupont White Teflon #807N-X; Magnet, ALNICO 5; Meet FDA Regulations 21CFR177.1550, 21CFR177.1520 and Teflon used to totally encapsulate the magnetic bar meets US Pharmacopeias Class VI. This product does not contain or come in contact Formaldehyde or other contaminants.
Sterilization	The products are non-sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	The product is Latex free.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	No special storage conditions are required when the goods are in transit. Suggested long term storage conditions are in an indoor, dry (laboratory) space with ambient temperatures of 65 C degrees (between 55 and 75C range), with low humidity (<50%), and away from excessive light and harsh chemicals.
Shelf Life	There is no definite shelf life or expiration for the product.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	February 2024
2	1.1	Updated Storage statement	May 2024