

Certificate of Quality

VWR Catalogue Number	77322-586
Description	VWR® Disposable Cleanroom Socks, White Color, Clothing size one, Unit Pack of 20 (1/ltn/Bx), Case of 400 (1/ltn/Bx)
Country of Origin	Manufactured in Korea
Date of Issue (yyyy-mm-dd)	2025-02-24

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with IEST-RP-CC003.4, Appendix B, Section B.2.5 (Helmke Drum Method).

Product Specifications

Material	Blend of Polyester / Nylon.
Sterilization	The product is non-sterile.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.

BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store at Room temperature.
Shelf Life	5 years

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	February 2025