

Certificate of Conformity

VWR Weighing Boat

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 611-9195

Lot No: 0840000000002

Capacity (ml)	Length (mm)	Pk
6	85	3

Material meet the following Specification

Acc. Internal Specification

Borosilicate Glass Glass 3.3 ISO 3585: Water resistance tested to ISO 719 (at 98 °C): HGB 1; Alkali resistance tested to ISO 695: Class A2; Acid resistance tested to ISO 1776: Na₂O<100 µg per 1 dm²

Quality Management DIN ISO 9001

This document has been produced electronically and is valid without a signature

Date of issue: **April 2020**