

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Critical Swabs
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2024-10-16

Catalogue Number	Description	Shaft/Head	Handle length [mm]	Pk
149-0330	Absorbent autoclavable swabs	PS/cotton	152	12
149-0331	Absorbent autoclavable swabs	Wood/cotton	152	1.000
149-0332	Large absorbent, autoclavable swabs	Wood/cotton	152	500
149-0340	Double-headed, cone-shaped, highly absorbent swabs	Paper/cotton	82	1.250
149-0341	Double-headed, highly absorbent swabs	Paper/cotton	78	2.500
149-0342	Double-headed, highly absorbent slim swabs	Paper/cotton	76	2.500

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 13485:2016**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Cotton tips, paper handles.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE)..
Cytotoxicity	No Information.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routine analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.
Storage Conditions	Dry.
Shelf Life	No expiry.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Annex

Specification:

Model	Handle Length [mm]	Handle Width [mm]	Head Length [mm]	Head Width [mm]	Total length [mm]
149-0330	152	2.5	16	5	168.0
149-0331	152	2.5	16	4.8	168.0
149-0332	152	2.5	18	6	170.0
149-0340	82	2.5	11	4.7	93.0
149-0341	78	2.5	15	4.4	93.0
149-0342	76	1.5	14	3.2	89.5

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	10-2025