

Material	330902Q
Material description	Ringer's solution tablets (quarter strength)
Grade	for milk testing

Characteristics	Specifications
Appearance	White tablets
Shape	Round biconvex
Clarity of solution	Passes test
Weight uniformity:	.
Max. 10 % outside 1150–1337 mg	Conforms
None outside 1120–1368 mg	Conforms
Average tablet weight	1181 - 1306 mg
Diameter	14.0 - 14.5 mm
Calcium (per tablet)	9.7 - 11.7 mg
Chloride (per tablet)	700 - 750 mg
Hardness	≥ 4 kp
pH	7 - 8

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.