



Product Specification Sheet

VWR SQUARE QUARTZ MICRO CELL, 100µl, LIGHT PATH 10 mm

VWR Cat.No.	VWRI634-0688
Product Name	SQUARE.QUARTZ MICRO CELL (10mm)
Serial Number	634-0688
Material	ES Quartz Glass
Material Thickness	1.25±0.1 (mm)
Size	12.5*12.5*45 L*W*H (mm)
Optical Path Length	10±0.05 (mm)
Optical Path Width	2±0.1 (mm)
Cubic Measure	100 (µl)
Transmittance	≥79% (at 200nm)
Package	1pcs/box
Beam Height	15mm

VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Date of issue: Dec.03rd 2012