



## Certificate of Conformity

---

# VWR ANTI-UV (2-1.2), ANTI-SCRATCH, ANTI-FOG, COVER SPECTACLES, LISBON

---

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

**Cat. No. :** 111-0930

**Type:** Clear PC

**Colour :** Transparent

**Weight:** 46 g

**Lens Thickness :** 2 mm

**Lens Base Curves :** 4

**Bridge Width :** 30 mm

**Bridge Height :** 20 mm

**Overall Width (hinge to hinge):** 150 mm

**Overall Length (lens- tips) :** 160 mm

### Other Characteristics:

- Visitor spectacles with wrap-around clear polycarbonate lens and with indirect ventilation.
- Anti-scratch/anti-fog coating
- Integrated brow guard
- Vented side shields
- Compliance to EN 166, EN 170; 1 FT norms

**Disclaimer:** VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2015**