



## ZPACK15150125 PLASTIC BOTTLE LIQUID 125 ML 24mm



### SPECIFICATIONS

|                   |          |
|-------------------|----------|
| Material :        | HDPE     |
| Neck Type :       | 24 mm    |
| Filling Opening : | 24,00 mm |
| Diameter :        | 1,986"   |
| Height :          | 3,868"   |
| Weight :          | 13,45 g  |

Remark : Packaging materials and specifications shown may be modified without preceding notice  
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011  
fax +32(0)16 385 385  
[info@be.vwr.com](mailto:info@be.vwr.com)  
[www.vwr.com](http://www.vwr.com)

Researchpark Haasrode 2020  
Geldenaaksebaan 464 B-3001 Leuven  
RPR/RPM Leuven  
BTW/TVA BE 0403.593.343