

Certificate of Conformity

VWR Catalogue Number	See List
Description	Magnetic stirring bars, Oval
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2026-06-24

Catalogue Number	Length [mm]	Diameter [mm]
442-0401	10	5
442-0402	15	6
442-0403	20	10
442-0404	25	12
442-0405	30	16
442-0406	35	16
442-0407	40	20
442-0408	50	20
442-0409	64	20
442-0410	70	20

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Virgin Polytetrafluoroethylene (PTFE) Grade TF1620.
Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI plastic (121°C).

Temperature Resistance -200 °C to +280 °C.

Sterilization	Non-Sterile.
Nitrosamine	Product is nitrosamine free.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	No special storage conditions required.
Shelf Life	Indefinite when stored under ideal conditions.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	06-2021
2	1.1	Additional information with regards to material has been updated	11-2022
3	1.2	Nitrosamine statement added	06-2023
3	2.0	Modified Quality system compliance and disclaimer	07-2025
4	2.1	Modified USP Class VI information	06-2026