

SECTION 1: Identification of the substance/mixture and of the company/undertaking

Identification of the product

Product name : Steroids Depleted Human Processed Serum
Reference : 33001D
Use of the substance/mixture : Raw material for research or manufacture of in vitro diagnostic reagents. Not for therapeutic use.

Identification of the company

Company : H2B
ESTER Technopole - 6 Allée Skylab
87068 LIMOGES Cedex- FRANCE
Contact : ☎ +33 5 19 08 00 20
☎ +33 5 55 35 95 44
✉ info.ivd@seqens.com

SECTION 2: Hazards identification

This substance/mixture is neither classified as dangerous for the Directive 67/548/EEC, nor classified in CLP Regulation n°1272/2008. Even if this product of human origin has been tested, it should be considered as potentially infectious and handled accordingly.

SECTION 3: Composition/information on ingredients

Nature of the product : Treated plasma, from human blood

Name	Product identifier	%	Classification
Human plasma	(CAS Number) - (EC Number) -	-	-

SECTION 4: First aid measures

Skin contact : Rinse out with water and soap.
Eyes contact : Rinse out with water. Seek medical advice.
Ingestion : Rinse out the mouth with water. Seek medical advice.
Inhalation : If irritation or sensitization, seek medical advice.

SECTION 5: Firefighting measures

Suitable fire extinguish media : Water + foam, CO2, powder
Non suitable media : None
Special exposure hazards : Do not breathe combustion products
Protection measures : Respiratory protective equipment

SECTION 6: Accidental release measures

Individual precaution : Wear protection gloves.
Environmental precautions : Do not disperse in the environment.
Cleaning procedures : Immediate cleaning with an absorbent matter (ex: Sawdust, sand) and clean up with water.

SECTION 7: Handling and storage

Handling : Handle with care to avoid microbial contamination.
Storage : Keep in the original packaging, between -16°C and -20°C and away from light.

Steroids Depleted Human Processed Serum

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

E007 IVDVP
Model B – June 2019
Version A – June 2019
Author : M.GAUDY

SECTION 8: Exposure controls/personal protection

Recommended personal protective equipment

Respiratory protection	: Respiratory protection not required.
Hands protection	: Wear protecting gloves (disposable).
Eyes protection	: Wear protecting glasses.
Body protection	: Impervious clothing. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.



SECTION 9: Physical and chemical properties

Appearance	: Limpid liquid
Colour	: Yellow to greenish yellow
pH	: 7.0 – 7.5
Solubility	: No data available
Boiling point	: No data available
Flash point	: No data available
Explosive properties	: No data available
Vapour pressure	: No data available
Division coefficient n-octanol/water	: No data available
Viscosity	: No data available
Vapour density	: No data available
Evaporation rate	: No data available

SECTION 10: Stability and reactivity

Stability	: Stable in recommended conditions of storage
Conditions to avoid	: Contaminated environment
Materials to avoid	: Strong oxidizing agent
Hazardous decomposition products	: Non available data

SECTION 11: Toxicological information

Acute toxicity	: No data available
LD/LC 50 data	: No data available
First irritating effect	
off the skin	: /
off the eyes	: /
Sensibilisation	: No sensitizing effect known
CMR effects	: Not classified CMR

SECTION 12: Ecological information

Ecotoxicity	: Non available data
Mobility	: Non available data
Degradability	: Non available data
Bioaccumulation potential	: Non available data
Others	: Do not disperse in the environment

Steroids Depleted Human Processed Serum

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH)

E007 IVDVP
Model B – June 2019
Version A – June 2019
Author : M.GAUDY

SECTION 13: Disposal considerations

Product	: Disposal according to the local administration advice, respecting laws and rules in force
Recommendations	: Eliminate the product in an approved centre of collection. Do not disperse in the environment.
Packaging	: Disposal according to the local administration advice. If there is no existing prescription thanks to a rule in force, dispose of as unused product.

SECTION 14: Transport information

Land Transport ADR / RID	: Not classified. Not hazardous product.
Sea Transport IMDG	: Not classified. Not hazardous product.
Air Transport OACI / IATA	: Not classified. Not hazardous product.

SECTION 15: Regulatory information

Labeling	: /
R-Phrases	: /
S-Phrases	: /

SECTION 16: Other information

This safety data sheet should be used in conjunction with technical sheet. The information given is based on our product knowledge, at the publication date of this document and is not exhaustive. The user must keep in mind the possible risks linked to different product uses than those for which it intended. It is the sole responsibility of the user to take all measures required in handling the product and, to know and apply the regulation concerning its activity. This does not alleviate the user from ensuring that legal obligation, different from those mentioned concerning use and storage of the product, do not exist.