

VWR International

Attn. Mr.Hadhi

Date : 22.10.13
Order-Nr. : e-mail 18.10.13 / Saleema Begam Abdul
Hadhi

Analysis-Certificate

Spoon spatula, Art. No.: 021.210


Typ of Material: 1.4301

C %	SI %	Mn %	P %	S %	Cr %	Ni %	
0,042	0,490	1,230	0,030	0,001	18,170	8,060	

The results only represent the result if one single analysis. These values vary with the individual material orders in the context of the limits put down on key to steel.

Best regards,

RSG Rostfrei-Schneidwerkzeuge GmbH



S.Mügge-Dewey



i.A. Simone Popko