

Quality Assurance Certificate

Arium® Scientific Kit

Order no. VWR-S-PACK
Use before 04 / 2027
Lot no. 2511483

This product is developed, produced and distributed according to a Quality and Environment Management System that is certified for compliance with ISO 9001 and ISO 14001.

The product is manufactured to meet the published specifications.

Statements

Hazardous Material:

The product does not contain any hazardous material as supplied and can be disposed in accordance to the local regulations.

Full Traceability:

The full traceability of the product is given in relation to all relevant information on raw materials, processing and handling.

Water Compliance:

The semiconductor (i. E. Electronics) grade ion exchange resin supplied in this cartridge kit is certified to be new and to produce 18.2 megohms + cm ultrapure water at 25°C according to meet the current published ASTM D1193 standard for ultrapure water.

Batch testing

Housing Leakage Test:

Each device has been individually tested for tightness by means of pressure drop test. The test has been performed according to Sartorius Standard Operation Procedures using a Sartocheck® testing unit. This device has passed the leakage test performed at a test pressure of 3 bar / 44 psi.

Visual Inspection:

100% of these products are inspected and pass SARTORIUS specifications.

Product Conformity:

Technical drawing compliance and batch record review.

06.11.2025

Date



Patricia Spannagel
Plant Quality Manager



Sartorius Lab Instruments GmbH & Co. KG
37070 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius.com