

Certificate of Conformity

VWR Catalogue Number	See List
Description	Safety caps
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2021-10-08

Catalogue Number	Description	Thread	Connector
590-0036	Safety cap I, including 1 PFA fitting and air valve	GL 45	1x for tubing 2,3/3,2 mm OD
590-0037	Safety cap II, including 2 PFA fittings and air valve	GL 45	2x for tubing 2,3/3,2 mm OD
590-0038	Safety cap III, including 3 PFA fittings and air valve	GL 45	3x for tubing 2,3/3,2 mm OD

Air Valve

Catalogue Number	Description
590-0035	Air valve, UNF 1/4", service life 6 months

Fittings for Capillary Connector Perfluoroalkoxy (PFA) material

Catalogue Number	Description	Colour
590-0042	Fitting with integrated ferrule, 1,6 mm OD	Green
590-0041	Fitting with integrated ferrule, 2,3 mm OD	Grey
590-0040	Fitting with integrated ferrule, 3,2 mm OD	Black
590-0039	Blind plug for capillary connector	Colourless

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	PTFE (Polytetrafluoroethylene) or PP (Polypropylene) Cap PP (Polypropylene) Valve PFA (Perfluoroalkoxy) fitting
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.

Storage Conditions Not Applicable.

Shelf Life Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.