



TECHNICAL DATA SHEET

Biohazard collection bags are printed with biohazard warning symbol and precautions in English, Spanish, and French. The included indicator patch darkens when autoclaved. Disposable.

Product name:	Biohazard Disposal Bags
Part number:	HS10320
Color:	Red
Thickness:	1.57 mil
Size:	7.99 x 12.00 inches (20.3 x 30.5 cm)
Autoclavable:	Yes to 250°F (121°C)

Heathrow Scientific®, LLC

620 Lakeview Parkway • Vernon Hills • Illinois • 60061 • Phone: (847) 816-5070 • Fax: (847) 816-5072
E-mail: info@heathrowscientific.com • Website: www.heathrowscientific.com