


## Recombinant Mouse CD34 Protein (Fc Tag) (A317476)

### Specifications:

|                    |                                                                                                                                                                                        |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name:              | Recombinant Mouse CD34 Protein (Fc Tag)                                                                                                                                                |
| Applications:      | SDS-PAGE                                                                                                                                                                               |
| Expression System: | HEK293 cells                                                                                                                                                                           |
| Nature:            | Recombinant                                                                                                                                                                            |
| Protein Species:   | Mouse                                                                                                                                                                                  |
| Protein Length:    | Protein fragment.                                                                                                                                                                      |
| Sequence:          | CD34(Thr35-Thr287)+hFc(Glu99-Ala330)                                                                                                                                                   |
| Tag:               | C-terminal Human Fc Tag                                                                                                                                                                |
| Molecular Weight:  | The protein has a predicted molecular mass of 53.3 kDa after removal of the signal peptide. The apparent molecular mass of mCD34-hFc is approximately 70-100 kDa due to glycosylation. |
| Conjugate:         | Unconjugated                                                                                                                                                                           |
| Purity:            | > 95% by SDS-PAGE and Coomassie blue staining.                                                                                                                                         |
| Product Form:      | Lyophilized                                                                                                                                                                            |
| Concentration:     | Reconstitution dependent.                                                                                                                                                              |
| Formulation:       | Lyophilized from sterile Phosphate Buffered Saline, pH 7.4. Normally 5%-8% Trehalose is added as a protectant before lyophilization.                                                   |
| Storage:           | Shipped at 4°C. Lyophilized: Store at -20°C to -80°C. Reconstituted: Aliquot and store at -80°C. Product is stable for one year. Avoid freeze/thaw cycles.                             |
| Disclaimer:        | This product is for research use only. It is not intended for diagnostic or therapeutic use.                                                                                           |

### Images:



Recombinant Mouse CD34 Protein (Fc Tag) (A317476) on SDS-PAGE under reducing conditions.