



Item Number IS115
Item Isopropyl Alcohol, 99%, USP, EP, BP, JP
CAS Number 67-63-0
Molecular Formula C₃H₈O
Molecular Weight 60.10

Assay (C₃H₈O)..... Min. 99.0%
Specific Gravity @ 25°C..... 0.783-0.787
Specific Gravity @ 20°C..... 0.785-0.788
Refractive Index @ 20°C..... 1.376-1.378
Acidity..... To pass test
Acidity or Alkalinity..... To pass test
Appearance..... To pass test
Clarity of Solution..... To pass test
Absorbance..... To pass test
Distilling Range..... 81°-83°C
(Min. 94% volume)

MAXIMUM LIMITS

Nonvolatile Residue..... 20 ppm
Benzene and Related Substances..... To pass test
Peroxides..... To pass test
Water..... 0.5%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843