

PRODUCT SPECIFICATION/CERTIFICATE OF CONFORMITY

PRODUCT CODE	1470/16D
DESCRIPTION	CRYSTALLISING DISH, FLAT BOTTOM
MANUFACTURER	PYREX®
COMPLIANCE	COMPLIES WITH DIN12337
MATERIAL	BOROSILICATE GLASS TYPE 3.3 WHICH COMPLIES WITH ISO3585 EUROPEAN AND US PHARMACOPOEIA GLASS
CAPACITY	2000ML
DIAMETER	190MM
HEIGHT	90MM
PACK QUANTITY¹	10

This is to certify that the goods detailed above have been inspected and tested in accordance with our BS EN ISO9001 and BS EN ISO13485 quality systems, and conform in all respects to the acceptable quality levels (AQL) detailed therein and any relevant specifications or drawings, unless recorded below as concessionary goods.

Signed: P Burton
Quality Technician
Date: 24 June 2021