



CERTIFICATE OF ANALYSIS

Product Description: Ansell-Kimtech Sterling Nitrile Gloves 24 cm Ambidextrous

Catalog Numbers: 99210, 99211, 99212, 99213, 99214

Issued Date: 8-Aug-2025

Total Cases per Lot #: 25	Lot: 2507UY	Batches: 25070907UY	Date of Manufacture: 2025-07-09	Expiration date: 2030-07-08
Total Cases per Lot #: 51	Lot: 2507UY	Batches: 25071207UY	Date of Manufacture: 2025-07-12	Expiration date: 2030-07-11
Total Cases per Lot #: 320	Lot: 2507UY	Batches: 25071755UY	Date of Manufacture: 2025-07-17	Expiration date: 2030-07-16
Total Cases per Lot #: 329	Lot: 2507UY	Batches: 25071753UY	Date of Manufacture: 2025-07-17	Expiration date: 2030-07-16
Total Cases per Lot #: 135	Lot: 2507UY	Batches: 25071852UY	Date of Manufacture: 2025-07-18	Expiration date: 2030-07-17
Total Cases per Lot #: 196	Lot: 2507UY	Batches: 25071855UY	Date of Manufacture: 2025-07-18	Expiration date: 2030-07-17

Inspection Test Data

	Watertight	Dimensions	Visual Defects		Physical Properties (Before Aging)	
			Minor	Major	ASTM - Ultimate Elongation (%)	ASTM - Tensile Strength (MPa)
Sample Size:	765	120	650	765	80	80
AQL Level:	0.65	4.0	4.0	2.5	2.5	2.5
Failures Allowed per AQL:	7	6	36	28	4	4
Failures:	0	0	0	0	0	0
Inspection Results:	Accept	Accept	Accept	Accept	Accept	Accept

Test Methods: Watertight (ASTM D 5151/EN 455-2), Elongation (ASTM D 412/ ISO 11193)/Tensile ASTM D 412/ FAB (EN 455-2/ISO 11193), ASTM D 6319

Physical Properties and Residual Powder Test

	Physical Properties (Before Aging)		Residual Powder
	ASTM - Ultimate Elongation (%)	ASTM - Tensile Strength (MPa)	(mg/glove)
Average:	555	39.72	0.04
Median:	551	39.48	0.04
Standard Deviation:	31	2.31	N/A
Minimum Limit:	500	28	2.00
Maximum Limit:	N/A	N/A	N/A
Inspection Result:	Accept	Accept	Accept

Test Methods: Elongation (ASTM D 412/ ISO 11193)/Tensile ASTM D 412/ FAB (EN 455-2/ISO 11193), Residual powder (ASTM D 6124 & EN 455-3)

Reviewed & Approved By:

Name: Florie Chin
Title: Senior Specialist, QA
Date: 8-Aug-2025

This document is intended for customer use only. Any disclosure, copying, distribution, or other use of this document is strictly prohibited.