



Certificate of Conformity

VWR LABELS, SQUARE

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

VWR Cat.Number	Description	Dimension (mm)	Pk
817-5000	Labels, white strips for 1.5 ml tubes, suitable for laser printers	33×13	1.700
817-5001	Labels, white strips for 1.5 - 2.0 ml tubes	33×13	1.000
817-5003	Labels, white strips for racks and boxes, suitable for laser printers	67×25	600
817-5005	Labels, multi-coloured strips for 1.5 ml tubes, suitable for laser printers	24×13	2.125
817-5006	Labels, white strips for 0.5 ml tubes, suitable for laser printers	24×13	2.380
817-5007	Labels, white strips for 0.2 ml PCR tubes	20×6	1.000
817-5008	Labels, yellow strips for 1.5 - 2.0 ml tubes	33×13	1.000
817-5011	Labels, blue strips for 1.5 - 2.0 ml tubes	33×13	1.000
817-5013	Labels, red strips for 1.5 - 2.0 ml tubes	33×13	1.000
817-5014	Labels, white strips for 96/384 well plates, suitable for laser printers	38×6	3.120
817-5015	Labels, green strips for 1.5 - 2.0 ml tubes	33×13	1.000
817-5016	Labels, purple strips for 1.5 - 2.0 ml tubes	33×13	1.000
817-5019	Labels, multi-coloured strips for 96/384 well plates, suitable for laser printers	38×6	3.900

Other Characteristics:

- The labels and adhesive have been designed to withstand temperatures from -196 to +121 °C.
- No need to laminate.
- Use on tubes, racks and boxes.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2016**