


CATALOG NO.		3410-42	
SIZE NOM VOL. (L).....		10	
NOMINAL BRIM CAP. (L)....		13.6	
MATERIALS	BOTTLE.....	PC	
	CLOSURE.....	WHITE PP	
	CLOSURE LINER..	SILICONE RUBBER	
NECK FINISH.....		48mm DOUBLE START THR'D	
THREAD (T.P.I.).....		3	
NOM. WEIGHT W/ CLOSURE (gms).....		755	
		INCH	CM (REF)
A (NECK I.D.).....		1.49±.05	3.8±.13
B (HT. W/O CLOSURE).....		13.17±.2	33.5±.51
C (HT. W/ CLOSURE).....		13.26±.25	33.7±.64
D (SQUARE).....		10.05±.25	25.5±.64
W.T. (NOM. WALL THK).....		0.08	0.2

2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

THIS PRODUCT SPECIFICATION
SHEET CONTAINS CONFIDENTIAL
AND/OR PROPRIETARY INFORMATION
AND IS TO BE USED ONLY FOR THE
EXPRESSED PURPOSES AGREED UPON
BY AN AUTHORIZED OFFICER OF
THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

MADE	PMP	APP'D	PSS 3410-42	REV.
CK'D	JC	PWC 12/16		1216
DESCRIPTION:				

BIOTAINER, CARBOY, PC, LAB, 10L