

1000 µL Conductive Tips

Specification Sheet (Part/REF # 235904)



1000 µL Conductive Non-Sterile Non-Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS Microlab Prep
Tip Type	Standard
Tip Volume	1000 µL
Maximum Tip Capacity	1,000 µL
Filter	Non-Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	White
Tip Material	Polypropylene - PP
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Rack Material	Polypropylene - PP
Tip Recycling Code	5
Rack Recycling Code	5
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack

Tip Length	95.1 mm
------------	---------

Orifice Size	0.7 mm
--------------	--------

Case Size	Case of 3,840 tips (Blister 5 x 96 tips per rack)
-----------	---

Packaging Option	Blister 96 Tip
------------------	----------------

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/automated-liquid-handling/disposable-tips/1000-µl-conductive-tips>