



Product Specification Sheet

VWR COLLECTION FUNNELS, LONG STEM, BOROSILICATE 3.3 GLASS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Cat. N°	Description	Stem diameter	Stem length
221-0172	Funnels, long stem, rim diameter 65 mm	9 mm	150 mm
221-0173	Funnels, long stem, rim diameter 80 mm	9 mm	150 mm
221-0174	Funnels, long stem, rim diameter 110 mm	9 mm	180 mm

Intended uses: Laboratory applications

Material: Borosilicate 3.3 glass

Shape: Conical funnels with long stem

Compliant to: DIN EN ISO 4800, Borosilicate glass 3.3 - properties

Chemical Properties: Glass has been tested by the Central Glass / Ceramic Research Institute, Kolkata, India and the report is found compliant to characteristics of Boro 3.3 Glass which also complies to requirements of EN ISO 3585, art.4.1,4.2,4.3a, 4 and ASTM E438-92.

Glass classification: Hydrolytic class HGB 1 according to ISO 719
Hydrolytic class HGA 1 according to ISO 720

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **October 2013**