

**Corning Incorporated
Life Sciences**

Registered
ISO 9001:2008

Product Description

Catalog Number: 3365

Product Description: Corning® 96-well polypropylene, round bottom plate, without lid

Component Materials:
Plate - Virgin Polypropylene, meets *USP, Class VI* requirements for plastic containers and closures.

Product Dimensions:

Length of Plate	-	5.035 in.	Diameter of Well @ top	-	.270 in.
Width of Plate	-	3.365 in.	Height	-	.560 in.
Depth to the bottom of well	-	.445 in.	Working Volume per well with storage mat sealed	-	200 µL
Tolerances of Dimensions	-	+/- .010 in.			

Surface Characterization:

Surface is characterized to be hydrophobic and non-charged for solution based assay.

RNase/DNase Testing:

This product has been tested and is free of any detectable RNase/DNase contamination.

Bovine Spongiform Encephalopathy and Transmissible Spongiform Encephalopathy:

This product complies with the latest revision of EMEA/410/01 "Note for Guidance on minimizing the risk of transmitting animal spongiform encephalopathy agents via human veterinary medicine products" by virtue of all bovine derived material having been processed per specific conditions of section 6.4 of EMEA/410/01.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for accurate labeling, and correct product configuration.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 12