

Corning[®] Aspirator

A convenient, user-friendly device that facilitates the aspiration of liquids from lab containers

The Corning Aspirator is an aspiration tool used for the removal of liquids from a variety of lab vessels. The device allows precise fingertip aspiration control, and is a sensible alternative to the use of glass pipets for supernatant disposal.

Designed for a comfortable fit in either hand, the aspirator enables easy handling of liquids from a variety of collection devices such as tubes, dishes, bottles, T-flasks, and multiwell plates.

- ▶ Vacuum regulation controlled by a touch-sensitive button allows for quick and comfortable aspiration of liquids
- ▶ Ergonomic design makes the device perfect for long work sessions or repetitive aspiration tasks
- ▶ The aspirator comes with a rubber adapter for accessory attachment, as well as a single-channel adapter for 1 to 200 µL disposable pipet tips
- ▶ An 8-channel adapter is available for working with 96-well microplates

- ▶ Anti-drip technology prevents liquids leaking from the pipet tips, reducing contamination risks and the likelihood of collecting erroneous results
- ▶ Entirely autoclavable for maximum safety when working with hazardous materials

Ordering Information

Products may not be available in all markets.

VWR Cat. No.	Corning Cat. No.	Description	Qty/Cs
89091-050	4930	Corning Aspirator	1

Accessory (sold separately)


VWR Cat. No.	Corning Cat. No.	Description	Qty/Cs
29443-002	4931	8-channel adapter for the Corning Aspirator	1



8-channel adapter

Related Products

Products may not be available in all markets.



VWR Cat. No.	Corning Cat. No.	Description	Qty/Pk	Qty/Cs
29443-352	4860	Corning® 1 - 200 µL universal fit racked pipet tips, yellow, polypropylene, sterile	96	960
89136-127	4121	Corning DeckWorks 1 - 200 µL pipet tips, graduated, hinged racks, natural, polypropylene, sterile	960	3,840
89136-173	4144	Corning DeckWorks 1 - 200 µL pipet tip station, graduated, natural, polypropylene, sterile	480	4,800
29442-416	4020	Stripette® 2 mL serological pipets, polystyrene, bulk packed, sterile	50	1,000
29443-043	4486	Stripette 2 mL serological pipets, polystyrene, individually paper/plastic wrapped, sterile	50	1,000
29442-418	4021	Stripette 2 mL serological pipets, polystyrene, individually plastic wrapped, sterile	100	1,000
29442-462	9186	Costar® 2 mL aspirating pipets, polystyrene, without graduations, individually wrapped, sterile	1	1,000

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.

avantor[™]
delivered by **vwr**[™]

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by country, federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2023 Avantor, Inc. All rights reserved.