

Data sheet

article: measuring cylinder 5000 ml, class B, amber stain graduation
article no.: 220 01 97

material:

Schott Duran (borosilicate glass 3.3) in accordance with ISO 3585

chemical properties:

- water resistance for 98°C: class 1 in accordance with ISO 719
- water resistance for 121°C: class 1 in accordance with ISO 720
- acid resistance: class 1 in accordance with DIN 12116
- alkali resistance: class 2 in accordance with ISO 695

physical properties:

- maximum service temperature: max. 500°C
- transformation temperature: 525°C
- density: 2,23 g/cm³
- linear expansion coefficient: $3,3 \times 10^{-6} \text{ K}^{-1}$

nominal volume:

5000 ml

subdivision:

50,0 ml

tolerance:

± 37,5 ml

norm:

ISO 4788

Hirschmann Laborgeräte
GmbH & Co.KG

i.A.S. Rössle

i.A. S. Rössle
Quality Manager

Eberstadt, 11.08.2021