
Date of Certificate: 10/3/2025

Catalog Number: 36902-0000
Description: STERILEWARE,SCOOP,PS,2OZ
Lot Number: J397502
Date of Manufacture: 8/21/2025
Country of Origin: US

Material: High Impact Polystyrene (HIPS)
We confirm that the material manufactured with HIPS conforms to:
21 CFR 177.1640
21 CFR 178.3297

Gamma Irradiation: Irradiation exposure dose (avg.) = 23.5 min - 50.0 max kGy
ANSI/AAMI/ISO 11137 (latest revision), Method 1
No testing is done for Pyrogens / Endotoxins.

CoSP Report Number: 3822788

Product Shelf Life: N/A - Product should remain in its irradiated state as long as outer packaging remains completely intact.

BSE/TSE Statement: Based on the information provided to us, the HIPS resin used in the manufacture of the above product is BSE/TSE free.
During manufacture, this product does not come into contact with any products or processes that may pose a BSE/TSE risk.

Other:

Bel-Art Products hereby certifies that the above listed product has been manufactured and inspected in accordance with documented procedures and specifications, and that it conforms to all applicable requirements.

Dejia Savage
Dejia Savage
Quality Inspector

Certificate of Processing

STERIGENICS 75 Tilbury Road Salem NJ 08079
 TEL 856 935-1400 FAX 856 935-0542 www.sterigenics.com

R55480102

09/03/2025 12:48:46 GMT
 Page 1 of 1

Customer Name:	Bel-Art Products-NJ	Processing Facility:	Salem	Work Order #	3822788
P.O.#	MD77197			Sales Order #	3881908
	23.5-50.0 kGy	BelArt, Gamma Treatment	Irradiation Date/Time:	08/31/2025 17:09:00 GMT	

SO Line #	Qty	UOM	Customer Item Number	Customer Item Description	Customer Lot Number	Customer Load Number
101.000	240	CA	369020000(20)	STERILEWARE,SCOOP,PS,2OZ	J397502	NA
	Dose Map		125_10885	0C3 - 0E3 x 1.14		
102.000	210	CA	369300000(30)	STERILEWARE,SAMPLING SPATULA,P	J398351	NA
	Dose Map		125_11119	0C3 - 0E3x1.095		
103.000	100	CA	369040006(20)	STERILEWARE®,METAL DETECT SCOOP	J397014	NA
	Dose Map		125_10885	0C3 - 0E3 x 1.14		
104.000	30	CA	369300000(30)	STERILEWARE,SAMPLING SPATULA,P	J397895	NA
	Dose Map		125_11119	0C3 - 0E3x1.095		
105.000	128	CA	379440000	STERILEWARE MINI-TONGS,PS,BX/2	J396168	NA
	Dose Map		125_10774	0C3-TE6		
	708	CA	Total			

Quality Test Summary

Op#	Quality Test Description	Minimum Spec	Maximum Spec	Result	Pass/Fail	-----Signed By----- User	Date /Time
450.00	Minimum Dose	23.5 kGy	50.0 kGy	25.7 KGY	Pass	ECEPEDA	09/01/2025 03:20:49 GMT
		Reason Code Test				EMANUEL CEPEDA	
450.00	Maximum Dose	23.5 kGy	50.0 kGy	40.4 KGY	Pass	ECEPEDA	09/01/2025 03:21:07 GMT
		Reason Code Test				EMANUEL CEPEDA	

Sterigenics certifies that the materials listed above (as described by the Manufacturer) received the indicated doses within the precision and accuracy of the dosimetry system employed.

Electronically Signed By: LESLIE GLESSNER
 Reason: Work Order Completions

Date: 09/03/2025 12:48:20 GMT