



Product Specification Sheet

VWR MAGNETIC STIRRING BARS, PLAIN

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Length (mm)	Ø (mm)	Pk
442-0262	10	6	10
442-0263	13	8	10
442-0264	15	4.5	10
442-0265	20	6	10
442-0266	20	7	10
442-0267	25	6	10
442-0268	25	10	10
442-0269	30	7	10
442-0270	35	10	10
442-0271	40	8	5
442-0272	50	10	5
442-0273	60	7	5
442-0274	65	13	5
442-0275	70	10	5
442-0276	75	13	5
442-0277	80	10	5
442-0278	57	27	5
442-0279	108	27	5
442-0280	159	27	5

Material: PTFE covered, strong Alnico V magnetic core

Properties: Good chemical resistance



Product Specification Sheet

Temperature Resistance: -200° up to +250°C

Compliant to: FDA and USP VI

PTFE used in the manufacture these products are pure, virgin material which conforms to FDA and USP VI requirements and the following specified substance content. It does not contain, or use in its manufacture any Class 1 ODC according to the Clean Air Act 1990, PFOA, POP, DEHP and BPA, no animal material is used in the manufacture or is present in this material.

	Description	Content (%)
Group A	Cadmium and Cadmium compounds	Not detected
	Hexavalent Chromium compounds	Not detected
	Lead and Lead compounds	Not detected
	Mercury and Mercury compounds	Not detected
	Tributyl Tin Oxide (TBTO)	Not detected
	Tributyl Tins & Triphenyl Tins	Not detected
Group B	Polybrominated Biphenyls (PBBs)	Not detected
	Polybrominated Diphenyl ethers (PBDE)	Not detected
	Polychlorinated Biphenyls (PCBs)	Not detected
	Polychloronaphthalenes	Not detected
	Short Chain Chlorinated Paraffins	Not detected
Group C	Asbestos	Not detected
	Azo Colorants	Not detected
	Ozone Depleting Substances	Not detected
	Radioactive Substances	Not detected

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **April 2013**