

Material	24667.290
Material description	Oil of castor
Grade	GPR RECTAPUR
CAS Number	8001-79-4

Characteristics	Specifications
Solubility in methanol	Passes test
Acid value	≤ 1.0
Density (20/4)	0.950 - 0.970
Iodine value	80 - 95
Saponification value	170 - 190
Non-saponifiable matter	$\leq 0.8 \%$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.