

Material	21199.361
Material description	Ammonium acetate
Grade	AnalaR NORMAPUR for soil analysis
CAS Number	631-61-8
Molecular formula	H ₃ CCOONH ₄
Molecular mass	77.08

Characteristics	Specifications
Assay	≥ 98.0 %
pH (20°C; 5 %)	5.0 - 7.4
Heavy metals (as Pb)	≤ 2 ppm
Ignition residue (SO ₄)	≤ 50 ppm
Cl (Chloride)	≤ 10 ppm
NO ₃ (Nitrate)	≤ 10 ppm
SO ₄ (Sulphate)	≤ 10 ppm
Ca (Calcium)	≤ 5 ppm
Cu (Copper)	≤ 0.5 ppm
Fe (Iron)	≤ 2 ppm
K (Potassium)	≤ 1 ppm
Mg (Magnesium)	≤ 1 ppm
Na (Sodium)	≤ 1 ppm
Pb (Lead)	≤ 0.5 ppm
Zn (Zinc)	≤ 0.5 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.