



# Product Specification Sheet

## VWR SQUARE QUARTZ MICRO CELL, 200µl, LIGHT PATH 10 mm

|                     |                                   |
|---------------------|-----------------------------------|
| VWR Cat.No.         | VWRI634-0689                      |
| Product Name        | SQUARE.QUARTZ MICRO CELL ( 10mm ) |
| Serial Number       | 634-0689                          |
| Material            | ES Quartz Glass                   |
| Material Thickness  | 1.25±0.1 (mm)                     |
| Size                | 12.5*12.5*45 L*W*H ( mm )         |
| Optical Path Length | 10±0.05 (mm)                      |
| Optical Path Width  | 4±0.1 (mm)                        |
| Cubic Measure       | 200 (µl)                          |
| Transmittance       | ≥79% (at 200nm)                   |
| Package             | 1pcs/box                          |
| Beam Height         | 15mm                              |

VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Date of issue: Dec.03rd 2012