

Specification

Buffer Solution pH 6.865

277124

Physical Description:	Liquid
Product Code:	277124
Product Name:	Buffer Solution pH 6.865
Short Description:	usage: These buffer solutions are useful for the calibration of high precision laboratory pH meters
Headline Comment:	* standardized, according to DIN19266
Specifications:	<p>pH (30°C): 6.853 ± 0.01</p> <p>pH (20°C): 6.881 ± 0.01</p> <p>pH (25°C): 6.865 ± 0.01</p> <p>pH (10°C): 6.900 ± 0.01</p>
WGK:	nwg
Storage:	RT
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 0072250192 • Register court Darmstadt HRB Nr. 7340