

Certificate of Conformity

VWR Dispenser Tips Classic 1.25 ml, Individual Blister

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 613-1002

Lot No: 9669021/1

Properties: Human DNA, RNase, and ATP Free

Human DNA Free: Samples will show no evidence of human DNA with a detection level < 0.5 pg/µl.

RNase Free: Samples will show no evidence of RNase with a detection level of < 1x10⁻⁹ Kunitz.

ATP Free: Samples will show no evidence of ATP with a detection level of <10⁻¹² mmol/ µl.

Production Date: November 2020

Expiry Date: November 2025

This document has been produced electronically and is valid without a signature

Date of Issue: **March 2021**