



Certificate of Conformity

VWR CARDBOARD CRYOBBOXES AND DIVIDERS

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	WxDxH ext. (mm)	Compartments	Height (mm)	Pk
479-1437	133x133x32	81 (9x9)	25	1
479-1438	133x133x50	81 (9x9)	30	1
479-1439	133x133x50	100 (10x10)	30	1
479-1440	133x133x75	81 (9x9)	40	1
479-1441	133x133x100	81 (9x9)	65	1
479-1442	133x133x130	81 (9x9)	62	1

Colour: White

Material: Cardboard



Other Characteristics:

- Resistance to temperature -140 °C.
- Waterproof coating.
- For storage of cryovials at low temperatures.





Certificate of Conformity

479-1437





479-1438







Certificate of Conformity

479-1439





479-1440





Certificate of Conformity


479-1441





Certificate of Conformity

479-1442



Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **September 2016**