



Sample organization for the busy lab

Droplet® Sample Storage Tray

Lightweight and sturdy to hold up to 24 full specimen containers (specimen containers not included).

- Removable insert to allow outer shell to be used for storage or transportation of laboratory products
- Stackable with other trays to support as many as four high
- Integrated label holder and cards to enable easily changed identification
- Able to be cleaned with common disinfectants and germicides

Part No.	Wells	Well Dia.	L x W x H	Pkg.	Color
HS120172	24	58 mm	18.9 x 13 x 3.3 in (48 x 33.1 x 8.4 cm)	1/ea	White ☐




Use for storage
or transport of
laboratory products





Trays are supplied "flat-packed"
reducing the amount of space needed for storage



Use for storage or transport of laboratory products



Removable inserts



droplet™

Corrugated sample storage tray

Disposable sample organization and storage

Corrugated storage trays can be used with the included dividers for specimen containers or remove them for general lab product storage or transport.

- Made from corrugated polypropylene
- Moisture resistant, significantly reducing the chance of mold
- Lightweight and sturdy to hold up to 24 full specimen containers
- Stackable with other trays
- Removable dividers allow outer shell to be used for storage or transport of laboratory products
- Pack includes: 10 trays and 10 dividers

- Disposable

Part No.	Wells	Well Dia.	L x W x H	Pkg.	Color
HS120375	24	58 mm	16 x 10.24 x 3.5 in (40.6 x 26 x 8.9 cm)	10/pk	White ☐