

Certificate of Conformity

VWR Catalogue Number	See List
Description	Mob caps, PP
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-12-23

Catalogue Number	Colour	Ø [mm]	Pk
113-9195	Blue	520	1.000
113-9194	White	520	1.000
113-9196	Green	520	1.000
113-9197	Red	520	1.000
113-9203	Yellow	580	1.000
113-9204	Lilac	580	1.000
113-9199	White	580	1.000
113-9200	Blue	580	1.000

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**.

QC Testing

Representative products samples are collected and inspected according to Level S-4 sampling plan7per MIL-STD-105E.

	Test Method	MD	CD
Break Tensile Properties [N]	ISO13934-1	10 N	7 N
Elongation rate [%]	ISO13934-1	100%	100%

Product Specifications

Material	Polypropylene.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted and have been shown to be non-toxic.
Non-Pyrogenic Statement	Product is pyrogen free.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols not present.
DEHP Statement	DEHP not present.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routine analysis is performed.
REACH Statement	Products meet the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006.
Storage Conditions	Room temperature, dry condition and shade from sun light.
Shelf Life	5 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	12-2025