



ACROS ORGANICS part of Thermo Fisher Scientific

ACROS
ORGANICS

Product Specifications

Product 41676
Papain, for biochemistry

General Product Data

Version	4
CAS No	9001-73-4
Molecular weight	0.00
Molecular formula	
Linear formula	

Product Specifications

Appearance (Color)	White to light yellow to light brown
Appearance (Form)	Powder
Infrared spectrum	Conforms
Assay	≥ 1000 TU/mg (Enz.act.) accuracy: $\pm 5\%$
Loss on drying	$\leq 7.0\%$ (60°C, 4 hrs) (vacuum)
pH	4.8 to 6.2 (1 in 50 solution)

Product Specifications are subject to amendment and may change over time.

Acros Organics

ENA23, zone 1, nr 1350, Janssen Pharmaceuticaalaa 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: <http://www.acros.com>
1 Reagent Lane, Fair Lawn, NJ 07410, USA
Fax 201-796-1329