

A10761 (1-Dodecyl)trimethylammonium bromide, 99%

Product Number: A10761
CAS number: 1119-94-4
MDL number: MFCD00011767
Molecular formula: C₁₅H₃₄NBr
Linear formula: CH₃(CH₂)₁₁N(CH₃)₃Br
Molecular weight: 308.35

Product Specification

Appearance (Color): White to cream
Form: Powder
Assay (Titration ex Bromide): >98.5 to <101.5%

Date of Print: December 10, 2020
Version: 7

Product Specifications are subject to amendment and may change over time.