

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Basic, Bouffant caps
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-07-23

VWR Catalogue Number	Diameter [mm]	Colour
113-8762	530	Blue
113-8879	580	Blue
113-8884	610	Blue
113-8874	530	Green
113-8880	580	Green
113-8885	610	Green
113-8877	530	Orange
113-8883	580	Orange
113-8888	610	Orange
113-8876	530	Red
113-8881	580	Red
113-8886	610	Red
113-1459	530	White
113-8878	580	White
113-1460	610	White
113-8875	530	Yellow
113-8882	580	Yellow
113-8887	610	Yellow

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 and ISO 13485:2016**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Inspected according to ISO 13934-1 standard.

	MD	CD
Tensile Strength [N/ 5 cm]	> 15	> 8
Elongation rate [%]	> 30	> 30

Product Specifications

Material	Spunbonded Polypropylene (SPP) 10g/m ² .
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.

RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	5 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.1	Created General CoQ using new layout	12-2021
2	2.0	Modified Quality system compliance, disclaimer and shelf life	07-2025