

Material	22544.292
Material description	Casein
CAS Number	9000-71-9

Characteristics	Specifications
Proteins (on anhydrous product)	≥ 90.0 %
Free acidity	≤ 1.0 %
Fats	≤ 2.5 %
Ignition residue (SO ₄)	≤ 2.75 %
Loss on drying	≤ 12 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.