



Product Specification Sheet

VWR TEST TUBES, POLYSTYRENE, NON-STERILE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Length (mm)	Ø ext.(mm)	Capacity (ml)	Pk
212-1837	75.0	11.8	3.5	5.000
212-1809	74.5	11.0	5	2.000
212-1834	74.5	12.8	5	4.000

Material: **Polystyrene, Clear**

Other Characteristics:

- Universal tubes with round bottom for a wide range of applications.
- Can be centrifuged up to 3000xg
- Can be delivered non sterile without plug
- Additional plugs and other colours upon request

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **November 2013**